



## HURLEY ADVANT 3/2MM WETSUIT CHEST ZIP

Code: MFS0013302-010

€ 174,96 ~~€ 249,95~~ - 30%

Hurley))

### HURLEY ADVANT 3/2MM WETSUIT

Hurley Advant 3/2 mm neoprene provides sustainable warmth and Exoflex stretch like no other. Crafted from premium Japanese limestone-based recycled neoprene for a premium feel and fit. Wetsuit liners and elastic webbing are dyed using a unique energy and water efficient process and strategically placed to add durability where you need it. Designed with insight from the world's top surfers and scientists at California State University San Marcos (USA) to ensure you can tackle colder waters without sacrificing comfort and flexibility.

Sustainable warmth: created using eco friendly processes including carbon black recycled from old tires, dope-dyed yarns, and limestone based neoprene.

Added stretch: New superstretch neospan exterior material allows you to move freely.

Quick dry, insulating interior: Xtend 2.0 interior lining provides you softer handfeel and faster drying.

Performance patterning: seams are placed out of your way so you can move freely.

Chest zip entry: watertight and easy to get on and off.

Key pocket placed at interior chest panel.

Spot taped : Double glued and blind stitched and strategically placed spot tape adds durability where you need it most.

Insulating interior: Back thermal wavy jacquard lining keeps you warm and cosy.

Smoothskin neck seal out water.

Lightweight kneepads: Kneepads feature jersey construction that absorbs less water.

Back right cuff Hurley Fastlane logo.

#### Surfcorner Store

Viale Umago 27 47921 Rimini, Italy

+39 0541 083880

Url: [https://www.surfcornerstore.com/hurley\\_advant\\_3\\_2mm\\_westsuit\\_chest\\_zip](https://www.surfcornerstore.com/hurley_advant_3_2mm_westsuit_chest_zip)



Filtri  
Size

**Surfcorner Store**

Viale Umago 27 47921 Rimini, Italy

+39 0541 083880

Url: [https://www.surfcornerstore.com/hurley\\_advant\\_3\\_2mm\\_westsuit\\_chest\\_zip](https://www.surfcornerstore.com/hurley_advant_3_2mm_westsuit_chest_zip)